



Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,

Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Sterile Eye Patch Specifications

Catalogue No. /Reference No. /Products Code	JC05025
Material	High-quality medical spunlace non-woven fabrics; Medical pressure-sensitive adhesives; Non stick absorbent pads; Release paper;
Color	White Color, Skin Color and Cartoon Printings
Size	5cm x 6.5cm, 5cm x 7.5cm, 5.8cm x 8.3cm, etc
Package	1 PC/Sterile Pouch
Shelf Life	5 Years
Sterile or Not	Sterile
Rules & Classification	Rule 4, Class I Sterile





Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,

Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Sterile Eye Patch Specifications

Reference No/ Product Code	Description	Size	Unit Package	Outside/Carton Package
JC050251	Sterile Eye Patch, Non Woven fabric, 1pc/sterile pouch, 20pouches/box, 50boxes/carton. EO Sterilization	7.5cm x 5cm White	1pc/pouch	20 pouches/box, 50 boxes/carton
JC050252		7.5cm x 5cm Skin Color		
JC050253		7.5cm x 5cm Cartoon Printing		
JC050254		5cm x 6.5cm White, Skin or Cartoon		
JC050255		5.8cm x 8.3cm White, Skin or Cartoon		